

950 Dimensions

	SIMPLICITY SEAT	DIMENSIONS
A	Minimum swivel radius With short armrest	26 ¾" 25 ¾"
в	Wall to inner seat back	3"
С	Seat depth	14 ½"
D	Minimum folded width to footrest	11 ¼"
E	Minimum open width to edge of footrest	22 7⁄8"
F	Armrest width - external	23 %"
G	Armrest width - internal	18 ¾"
H	Seat back height	15 ¾"
0	Footrest to seat height Standard seat option Extended seat post	15 ¾" – 19 ¼" 18 ⅛" – 22 ⅛"
J	Minimum footrest height	2"
	Minimum track intrusion into staircase	6 ¼"
	Maximum weight limit Manual hinge or Slide Track	310 lbs 265 lbs



Upholstery is made from a wipe clean and flame-retardant material.


All of our stairlifts are tested and approved by worldwide independent inspection services to give you the assurance that our products meet the highest standards in terms of quality and safety.







Handicare 10888 Metro Court St. Louis, MO 63043 USA 888-637-8155 sales.us@handicare.com www.handicareusa.com



STRAIGHT STAIRLIFT



## STRAIGHT STAIRLIFT

Our economical, entry level 950 stairlift is everything but. Not only does the 950 get you safely up and down the stairs, but it is designed to get you there with the highest level of convenience. Handicare has been a world leader in innovative lifting and mobility solutions for 130 years, providing safe solutions such as the 950 stairlift.

- Swivel your seat toward the landing for a safer exit when you reach the top
- Each unit comes with 2 handheld remotes to easily retrieve the stairlift
- Prevent unauthorized use with provided key
- Upholstery provided in a sand-colored wipe clean flame retardant fabric
- Limited lifetime warranty on the motor and gearbox; 2 year manufacturer's parts warranty

An innovative and economical solution for your straight staircase

# Ergonomic and innovative design features

### **SMART**



Stairlifts are battery powered to ensure you are in motion even during power outage. The continuous charging feature charges the battery no matter where it is parked.

**COMPACT** 

your seat.



With a foldable footrest, seat and armrests, the 950 boasts a slim profile of only 111/4 inches, ensuring other stair users can still use the staircase comfortably and safely.

Your first step to independence

### SAFE



Surrounded by safety sensors, Handicare stairlifts are designed to stop if meeting an obstruction - keeping you safe and still on

### **REVOLUTIONARY**



The Slide Track moves up the staircase at the same time the seat is moving, leaving the bottom of the stairs completely unobstructed. A manual hinge is also available.



### INTUITIVE



The toggle control was designed in collaboration with occupational therapists to ensure it is easy-to **use** and operate – simply push with your hand, arm or wrist in the direction you want to go.

### ADAPTABLE



The **zero intrusion** option is used when a doorway close to the top of the stairs would be obstructed by the track. (This option is not available with the hinge track option.)

